

## THC-U11

## Marijuana THC Rapid Test Strip (Urine)

### INTENDED USE

The THC Rapid Test Strip (Urine) is a rapid visual immunoassay for the qualitative presumptive detection of THC metabolites (11-nor- $\Delta^9$ -THC-9-carboxylic acid) in human urine specimens.

### INTRODUCTION

Marijuana, cannabis or tetra-hydro-cannabinol (THC) is a hallucinogenic agent derived from the flowering portion of the hemp plant. Smoking is the primary method of use of marijuana/cannabis. Higher doses used by abusers produce central nervous system effects, altered mood and sensory perceptions, loss of co-ordination, impaired short-term memory, anxiety, paranoia, depression, confusion, hallucinations and increased heart rate. A tolerance to the cardiac and psychotropic effects can occur, and withdrawal syndrome produces restlessness, insomnia, anorexia and nausea. When marijuana is ingested, the drug is metabolised by the liver. The primary urinary metabolite of marijuana is 11-nor- $\Delta^9$ -THC-9-carboxylic acid, and its glucuronide. This means that the presence of detected cannabinoids, including the primary carboxyl metabolite, in the urine indicate marijuana/cannabis use.

### PRINCIPLE

The THC Rapid Test Strip (Urine) has been designed to detect the THC metabolites through visual interpretation of color development in the strip. The membrane was immobilized with THC conjugates on the test region, and the sample pad was pre-coated with colored anti-THC metabolites antibodies colloidal gold conjugates. After specimens were added, the gold-conjugates move along the membrane chromatographically by capillary action and antibodies get to the test region. If there is no drug molecule in the urine the antibody gold conjugate would attach to the drug conjugate to form a visible line. Therefore, the formation of a visible precipitant in the test region occurs when the urine is negative for the drug. If THC metabolites are present in the urine, the drug antigen competes with the immobilized drug conjugate on the test region for limited antibody sites. In case of sufficient concentration of the drug, it fills the limited antibody binding sites. This will prevent attachment of the colored antibody-colloidal gold conjugate to the drug conjugate zone on the test region. Therefore, absence of the colored band on the test region indicates a positive result. Appearance of a colored band at the control region serves as a procedural control. This indicates that proper volume of specimen has been added and membrane wicking has occurred.

### REAGENTS

Each test consists of a reagent strip mounted in a plastic housing. The amount of each antigen and/or antibody coated on the strip is less than 0.001 mg for antigen conjugates and goat anti-rabbit IgG antibodies, and less than 0.0015 mg for antibody components. The control zone of each test contains goat anti-rabbit IgG antibody. The test zone of each test contains drug-bovine protein antigen conjugate, and the conjugate pad of each test contains monoclonal anti-drug antibody and rabbit antibody-colored particle complex.

### MATERIALS

#### Materials Provided

- Individually packed test strips
- Package insert

#### Materials Required but Not provided

- Specimen collection container
- Timer

### PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Do not use after expiration date indicated on the package. Do not use the test if its foil pouch is damaged. Do not reuse tests.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore, recommended that these products be treated as potentially infectious, and handled observing the usual safety precautions (do not ingest or inhale).
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to performing any tests.
- Do not eat, drink or smoke in the area where the specimens and kits are handled. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Humidity and temperature can adversely affect results.
- The used testing materials should be discarded in accordance with local, state and/or federal regulations.

### STORAGE AND STABILITY

- The kit should be stored at 2-30°C until the expiry date printed on the sealed pouch.
- The test must remain in the sealed pouch until use.
- Do not freeze.**
- Care should be taken to protect the components of the kit from contamination. Do not use if there is

evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to false results.

### SPECIMEN COLLECTION AND STORAGE

- The THC Rapid Test Strip (Urine) is intended only for use with human urine specimens.
- Collected urine specimens must be put in clear and dry containers. Ensure that a sufficient quantity of the specimen is collected to allow submerging the dipping area of the strip.
- Perform the testing immediately after the specimen collection. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 2-8°C for up to 48 hours. For long term storage, specimens should be kept below -20°C.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Avoid repeated freezing and thawing of specimens.
- Pack the specimens in compliance with applicable regulations for transportation of etiological agents, in case they need to be shipped.

### PROCEDURE

Bring tests, specimens, and/or controls to room temperature (15-30°C) before use.

- Remove the test from its sealed pouch and use it as soon as possible. To obtain a best result, the assay should be performed within one hour.
- Hold the strip at the handle with the product name imprints. Do not touch the membrane part of the strip to avoid contamination.
- Dip the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing the strip. As the test begins to work, you will see color move across the membrane.
- Take the strip out of the specimen afterwards and place it on a non-absorbent flat surface. Start the timer and wait for the colored line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.

### INTERPRETATION OF RESULTS



**POSITIVE: Only one colored band appears, in the control region (C).** No apparent colored band appears in the test region (T).



**NEGATIVE: Two colored bands appear on the membrane.** One band appears in the control region (C) and another band appears in the test region (T).



**INVALID: Control band fails to appear.** Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

#### NOTE:

- The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered negative. Note that this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.
- Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

### QUALITY CONTROL

- Internal procedural controls are included in the test. A colored band appearing in the control region (C) is considered an internal positive procedural control, confirming sufficient specimen volume and correct procedural technique.
- External controls are not supplied with this kit. It is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

### LIMITATIONS OF THE TEST

- The THC Rapid Test Strip (Urine) is for professional *in vitro* diagnostic use, and should be used for the qualitative detection of THC metabolites only.
- The THC Rapid Test Strip (Urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
- Please take the specificity and the cross reactivity into account for evaluation.
- A positive result with any of the tests indicates the presence of a drug/metabolite only, and does not indicate or measure intoxication.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- There is a possibility that technical or procedural errors as well as other substances and factors not listed may interfere with the test and cause false results.
- The test is designed for use with human urine only. Due to absence of ions and other components in pure water the usage of pure water for test could lead to false or invalid results.
- The test does not distinguish between drugs of abuse and certain medications.

### PERFORMANCE CHARACTERISTICS

#### A. Accuracy

The accuracy of the THC test was compared and checked against a commercially available test with a

cut-off value of 50 ng/ml. 120 urine samples taken from volunteer test persons who claim to be non-consumers was examined under both tests. The results were 100% in agreement.

#### B. Reproducibility

The reproducibility of the THC test was verified by blind tests performed at a four different locations. Of the 60 samples with 11-nor- $\Delta^9$ -THC-9-carboxylic acid concentration of 25 ng/ml, all were determined negatives. Of the 60 samples with 11-nor- $\Delta^9$ -THC-9-carboxylic acid concentration of 100 ng/ml, all were determined positives.

#### C. Precision

Test precision was determined by blind tests with control solutions. Controls with 11-nor- $\Delta^9$ -THC-9-carboxylic acid concentration of 25 ng/ml should yield a negative result and controls with 11-nor- $\Delta^9$ -THC-9-carboxylic acid concentration of 75 ng/ml should provide a positive result.

#### D. Specificity

The specificity of the THC Test was tested with the substances listed below, all of which can be found in a normal urine specimen.

The following compounds with a similar chemical structure yield a positive result at the specified concentration:

THC-test with Cut-off 50 ng/ml (THC50)

Drug	Concentration (ng/ml)
11-nor- $\Delta^9$ -THC-9-COOH	50
11-nor- $\Delta^9$ -THC-9-COOH	50
$\Delta^9$ -Tetrahydrocannabinol	15000
$\Delta^9$ -Tetrahydrocannabinol	20000
Cannabinol	20000

With exception of the above, for the respective parameter listed positive-reacting drugs resp. drug metabolites, all following listed compounds reacted negative up to a concentration of 100 µg/ml.

Acetaminophen	Guaiacol Glyceryl Ether
Acetone	Hemoglobin
Albumin	Imipramin
Amitriptylin	(+/-)-Isoproterenol
Ampicillin	Lidocaine
Aspartam	(+)-Naproxen
Aspirin	Oxalic Acid
Atropine	Penicillin-G
Benzocaine	Pheniramine
Bilirubin	Phenothiazine
Coffeine	Phenylethylamine
Chloroquin	Procaine
(+/-)-Chlorpheniramin	Quinidin
Chlorpheniramin	Ranitidin
Creatin	Riboflavine
Dexbrompheniramine	Sulindac
Dextromethorphan	Sodium Chloride
4-Dimethylaminoantipyrine	Thioridazin
Dopamine	Trifluoperazin
Erythromycin	Trimethobenzamid
Ethanol	Tyramine
Furosemide	Vitamin C
Glucose	

### LITERATURE REFERENCES

- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd ed. Davis: Biomedical Publications; 1982.
- Urine Testing for Drugs of Abuse. National Institute on Drug Abuse (NIDA), Research Monograph 73, 1986.
- Thomas L. eds., Labor und Diagnose, 6. ed., TH-Books Verlagsgesellschaft, Frankfurt, 2005.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53, 69, 11970, 1988.
- McBay AJ. Drug-analysis technology--pitfalls and problems of drug testing. Clin Chem. 1987 Oct; 33 (11 Suppl): 33B-40B.
- Gilman AG, Goodman LS, Gilman A, eds. Goodman and Gilman's The Pharmacological Basis of Therapeutics. 6th ed. New York: Macmillan; 1980

### GLOSSARY OF SYMBOLS

	Catalog number		Temperature limitation
	Consult instructions for use		Batch code
	<i>In vitro</i> diagnostic medical Strip		Use by
	Manufacturer		Do not reuse